


FastRunner DNA Ladder

(50bp - 2,000bp)

Catalogue Number: 12800

100 gel lanes





FastRunner DNA Ladder (50bp - 2,000bp)

Catalogue Number: 12800 (100 gel lanes)

High quality ready-to-use double-stranded DNA ladder derived by recombinant technology for precise on-gel size and mass determination.

- Precise band sizes and masses for densitometric quantification.
- Eight discrete bands, ranging from 50bp to 2,000bp.
- Higher intensity band at 500bp for easy reference.

1mL of premixed DNA ladder (0.5µg/10µL) in loading buffer (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, and 0.17% SDS).

Mix thoroughly. For best results, load 10µL of DNA ladder per well. For precise mass determination with a densitometer, stain gel after electrophoresis using 0.5µg/mL ethidium bromide for 30-40 minutes. The accompanying diagram shows the size and mass for each band based on 10µL ladder per well.

Stable at room temperature. For longer term storage, -20°C is recommended.

This ladder was standardized using 10µL of DNA per lane on a 0.8 cm thick, 13 x 15 cm, 2.0% agarose gel run in TAE buffer.

Tech support: techsupport@norgenbiotek.com

Orders: orders@norgenbiotek.com

Toll Free: 1 866 6674362

Fax: 905 227 1061

www.norgenbiotek.com