



Exceptional Transfection Efficiency in Human Induced Pluripotent Stem Cells (iPSCs) via Reverse Transfection with *TransIT*[®]-LT1. The *TransIT*[®]-LT1 Transfection Reagent was used to reverse transfect 1.3×10^6 iPS cells with a ZsGreen expressing plasmid (Clontech). Reverse transfections were performed in 6-well plates using 12 μ l of *TransIT*[®]-LT1 Transfection Reagent to deliver 4 μ g of DNA (3:1, reagent: DNA). Cells were visualized 48 hours post-transfection and imaged under a 10X objective with an Olympus IX71[®] Inverted Microscope. Images are (A) phase contrast and (B) green fluorescence. Cells were assayed 48 hours post-transfection on an Accuri[®] Cytometer. The histogram (C) shows untransfected cells (black line) compared to cells transfected with plasmid using *TransIT*[®]-LT1 (green line). [See more information on stem cell applications.](#)



High Efficiency Transfection of iCell[®] Cardiomyocytes Using *TransIT*[®]-LT1 Transfection Reagent. iCell[®] Cardiomyocytes were plated at 20,000 cells/well in a 96 well tissue culture plate coated with 0.1% gelatin. After allowing the cells to recover from thaw, cells were transfected with 100 ng/well of pMAXGFP (Lonza) using *TransIT*[®]-LT1 Transfection Reagent with a 2:1 reagent-to-DNA ratio according to the manufacturer's instructions. Fluorescent images were taken 3 days post transfection using a Olympus IX71[®] inverted microscope. [See more information on stem cell applications.](#)

Data courtesy of





The *TransIT*®-LT1 Reagent Exhibits Higher Expression and Lower Cellular Toxicity Compared to Other Transfection Reagents. HepG2 cells were transfected with a luciferase expression plasmid using the designated reagents at the manufacturer’s recommended reagent-to-DNA ratio indicated beneath each bar. Transfections were performed in 96-well plates using 0.1 µg of plasmid DNA per well. Luciferase expression (bar graph) and lactate dehydrogenase (LDH) levels (line graph) were measured at 24 hours post-transfection. LDH levels are reported as % cytotoxicity compared to cells alone and were measured using a commercially available colorimetric assay; all values at or below zero are represented as zero on graph. Experiments were performed as per industry accepted testing protocols. FuGENE® is a registered trademark of Fugent LLC. Lipofectamine™ is a trademark of Life Technologies Corporation.