


## Figures and Data

[Delivery of Cas9 mRNA and CRISPR Guide RNA with \*TransIT\*®-mRNA Transfection Kit High Efficiency and Low Toxicity Transfection Following 14 Consecutive Transfections with \*TransIT\*®-mRNA Transfection Kit](#)

[TransIT®-mRNA Transfection Kit Transfects GFP mRNA into DC 2.4 Dendritic Cells Multiple Dendritic Cell Types Express GFP from mRNA Transfected by \*TransIT\*®-mRNA Transfection Kit](#)


[High Level Luciferase Expression after Delivery of a Luciferase mRNA using \*TransIT\*®-mRNA Transfection Kit](#)

[TransIT®-mRNA Transfection Kit Efficiently Delivers the lacZ mRNA to CHO-K1 Cells](#)




### Efficient Genome Editing with Cas9 mRNA + Guide RNA Oligonucleotides.

HEK293T/17, U2OS and NHDF cells were co-transfected with 0.5  $\mu\text{g}$  of Cas9 encoding mRNA, 5meC,  $\psi$  (Trilink Biotechnologies) and 25nM of PPIB targeting 2-part gRNA (Dharmacon) using *TransIT*®-mRNA Transfection Kit (0.5  $\mu\text{l}$ /well of 24-well plate of both mRNA Reagent and Boost, Mirus Bio). A T7E1 mismatch detection assay was used to measure cleavage efficiency at 48 hours post-transfection.




**High Efficiency and Low Toxicity Transfection Following 14 Consecutive Transfections with *TransIT*®-mRNA Transfection Kit.** Repeated daily transfections were performed in the same population of BJ fibroblasts using *TransIT*®-mRNA Transfection Kit, Lipofectamine® RNAiMAX (Life Technologies) and Stemfect™ RNA Transfection Kit (Stemgent) - with a capped and polyadenylated EGFP mRNA incorporating pseudouridine and 5mC modified bases (Trilink Biotechnologies, Inc.). Multiple reagent-to-RNA ratios were tested and the optimal ratio is represented. Transfections were performed in 12-well plates using the indicated reagent-to-RNA ratios to deliver 1 µg of RNA. GFP and phase contrast images were taken in the same field of view everyday after transfection. Transfection efficiency was measured by flow cytometry on a Guava® easyCyte™ 5HT following 14 consecutive daily transfections (blue bars). Cell viability was determined using cell counts measured during flow cytometry (black line grey bars). Error bars represent the standard deviation of triplicate wells.



***TransIT®-mRNA Transfection Kit Transfects GFP mRNA into DC 2.4 Dendritic Cells.*** Using *TransIT®-mRNA* Transfection Kit, DC 2.4 cells were transfected with (A) 0.5 µg, (B) 1 µg and (C) 2.5 µg of capped and polyadenylated mRNA encoding GFP with 1 µl *TransIT®-mRNA* Reagent and 1 µl Boost. Cells were seeded overnight at 100,000 cells/well in 24-well plates. Images were taken 10 hours post-transfection.


*Data courtesy of Kyle Phua (Principal Investigator: Kam W. Leong), Duke University*

## Transfection Efficiency on Primary Murine BMDC




**Multiple Dendritic Cell Types Express GFP from mRNA Transfected by *TransIT*<sup>®</sup>-mRNA Transfection Kit.** Murine primary bone marrow derived dendritic cells (BMDC) and murine dendritic cells types (JAWSII and DC 2.4) were transfected with 1  $\mu$ g of capped and polyadenylated mRNA encoding GFP using a *TransIT*<sup>®</sup>-mRNA Reagent: Boost: mRNA ratio of 1:1:1 ( $\mu$ l: $\mu$ l: $\mu$ g). Primary BMDCs, JAWSII and DC 2.4 were seeded (80,000 cell/well) overnight in 24-well plates. Cells were assayed via flow cytometry 8 hours post transfection. Error bars represent the standard deviation of at least 3 separate experiments.


*Data courtesy of Kyle Phua (Principal Investigator: Kam W. Leong), Duke University*



**High Level Luciferase Expression after Delivery of a Luciferase mRNA Using *TransIT*®-mRNA Transfection Kit.** Cells in 12-well plates were transfected with a capped and polyadenylated mRNA encoding luciferase using the *TransIT*®-mRNA Transfection Kit. Approximately 18 hrs post-transfection the cells were harvested and the total luciferase activity per well was determined.



A. Mock Transfection Control



B. *LacZ* mRNA Transfection

***TransIT*®-mRNA Transfection Kit Efficiently Delivers *lacZ* mRNA to CHO-K1 Cells.** Using the *TransIT*®-mRNA Transfection Kit, CHO-K1 cells were mock transfected (A) or transfected with a capped and polyadenylated *lacZ* encoding mRNA (B). Approximately 18 hrs post-transfection the cells were stained using Mirus Bio's Beta-gal Staining Kit to identify the *lacZ* transfected cells.