



Figures and Data

[TransIT®-Oligo Reagent Achieves High Transfection Efficiency](#)

[TransIT®-Oligo Reagent Effectively Transfects a 2'OMe RNA Oligo that Blocks a Cryptic Splice Site](#)



TransIT®-Oligo Reagent Achieves High Transfection Efficiency. HeLa cells transfected using *TransIT®-Oligo Reagent* and *Label IT® Cy®3* and *Label IT® Fluorescein* labeled phosphothioate DNA oligos in complete media for 24 hours.



***TransIT®-Oligo* Reagent Effectively Transfects a 2'OMe RNA Oligo that Blocks a Cryptic Splice Site.** The HeLa-Luc 705 reporter cell line (Kang *et al.* 1998, 37:6235) used in this study contains a luciferase reporter construct that has the β -globin 705 intron inserted into the luciferase ORF. A mutation present at position 705 of the β -globin intron activates two cryptic splice sites within the intron that lead to the production of a spliced luciferase mRNA that is disrupted by a small intron with an in-frame stop codon, thus preventing translation of functional luciferase protein. The transfection of a 2'OMe oligonucleotide ([TriLink BioTechnologies, Inc.](#)) complementary to the cryptic 705 splice site inhibits splicing at the cryptic splice sites enabling the complete removal of the β -globin intron and production of a mRNA with a complete, uninterrupted luciferase ORF.

The HeLa-Luc 705 cell line was transfected with increasing amounts of the anti-705 splice site 2'OMe RNA oligo at the indicated final concentrations using the *TransIT®-Oligo* Transfection Reagent. The cells were harvested 24 hours post-transfection and assayed for luciferase activity. The increase in luciferase activity indicates effective delivery of the anti-705 splice site RNA oligo using the *TransIT®-Oligo* Reagent.