



***TransIT*[®]-3T3 Transfection Kit**

A high efficiency, low toxicity, DNA transfection reagent optimized for 3T3 cell types

- **Unique Formulation**—Each kit contains a base and booster reagent to maximize transfection performance in 3T3 cell types.
- **Low Cellular Toxicity**—Maintain cell density and reduce experimental biases.
- **High Efficiency Delivery**—Achieve 55-65% transfection efficiency in NIH 3T3 cells to ensure experimental success.

Data

Figure 1. The *TransIT*-3T3 Kit Achieves High Transfection Efficiency.

NIH 3T3 cells were transfected with pEGFP using the *TransIT*-3T3 Kit in complete media for 24 hours.

