

Kit Specifications	
Column Binding Capacity	50 µg
Maximum Column Loading Volume	650 µL
Time to Complete 10 Purifications	30 minutes
Maximum Amount of Starting Material: Fungi (wet weight) Yeast or Bacteria Culture	50 mg 1 x 10 ⁸ cfu's
Average Purity (OD260/280)	(~0.5 - 1 mL)
Average Yield* <i>Pichia sp.</i> (yeast) <i>Botrytis cinerea</i> <i>Fusarium sp.</i> <i>Aspergillus niger</i> <i>Penicillium sp.</i> <i>Alternaria tenuissima</i>	25 µg 32 µg 42 µg 26 µg 40 µg 30 µg

* Yield will vary depending on the type of sample processed