

Parameter	Gel Extraktion	PCR Clean-up
Max. sample volume	300 mg Agarosegel	100 µl PCR Mix
Gel	< 2,5% TAE oder TBE	
Typical recovery	70-80%	80-90%
Binding capacity	10 µg	10 µg
DNA fragment length	50 bp – 10 kb	50 bp – 10 kb
Primer removal		< 25 bp
Elution volume	20 – 50 µl	20 -50 µl
Prep time	20 minutes	20 minutes
Format	spin column	spin column