

Sample Hold (volume of sample per 10ml Ecoscint XR)											
Water										10ml	
Sodium Hydroxide 0.1N								8ml			
Sodium Chloride 0.5M					5ml						
Sucrose 32%											
Ammonium Acetate 15mM										10ml	
Ammonium Acetate 0.25M										10ml	
Tris HCl 0.05M pH 7											>10ml
Tris HCl 0.5M pH 7										10ml	
Washing Buffer 1X PBS										10ml	
Ammonium Sulfate 0.1M										10ml	
Methylene Chloride	All proportions										
Acetonitrile	All proportions										
Ammonium Formate 1M					4ml						