

Number of gels	1-4
Precast gel compatibility (Up to 2 gels/run)	IDGel, Novex, SERVAGel, Thermo Precise Pierce Protein Gel
Plate dimensions (w x h x t)	10x10x0.2cm
Gel Dimensions (w x h)	8x8.5cm
Total buffer Volume for 2 gels	Min: 250ml; Max: 1200ml
Total buffer Volume for 4 gels	Min: 250ml; Max: 1200ml
Standard run time for SDS-PAGE	1-2 hours at from 90-225V
Recommended power supplies	P5-0051; P5-0048; P5-0054for blotting
Unit Dimensions (w x d x h)	19x13x15cm
Weight	1.8Kg