

# Fermacidal D2

Disinfectant for surfaces,  
objects and tools

*For disinfection and contamination prevention of laboratory devices and furniture. Suitable for use on most surface materials, non-corrosive, not spotting.*

*Keeps your cultures safe - harmless for cells and tissues, non-volatile, non-toxic, non-irritant.*



**Ingredients:** Quarternary ammonium compounds. Does not contain alcohol, mercury, formaldehyde, phenol or other toxic compounds.

**Efficiency:** Bactericidal (*E. coli*), fungicidal (*Candida*), efficient against mycobacteria (*TbB*), mycoplasma, hepatitis B, HIV, and other viruses (e.g. rota, influenza A, H5N1, H1N1, corona, etc.), prevents endospores from flourishing.

**Properties:** Odourless, non-volatile, harmless for mammalian cells and tissues in culture, non-irritant, skin-friendly, biodegradable.

**Instructions for use:** Cell culture flasks and plates may be kept in the incubator during application. Spray objects daily until wet. Then let surfaces air-dry. Drying time = Reaction time, no wiping!






**Available sizes:** Available in sizes from 50 ml up to 10 L.

**Safety information:** Disposal of the biocidal product and its packaging according to official regulations. P102 – Keep out of reach of children. Harmful if not used as directed.

**Application:** For cleaning, disinfection and contamination prevention of surfaces, objects and tools (cell and tissue culture incubators, clean benches, furniture, devices, instruments, tools, tables, handles, handrails, tapware, etc.) in cell culture and molecularbiology laboratories, hospitals, therapy- and nursing facilities, etc.

**Registered as biocide:**

BAG / FOPH No. CHZB0445  
BauA Reg.-No. N 24040  
VAH-listed

Type	Spray			Container	Cloths
List					
Size	50 ml	250 ml	1000 ml	5 l / 10 l	120 pcs.
Part Code	215.001	217.001	219.001	222.005 (5 l) 222.010 (10 l)	225.001 (box) 226.001 (refill)

**Ordering information:**



**GeneFlow Ltd**  
1 The Sycamore Tree  
Elmhurst Business Park  
Lichfield  
Staffordshire  
WS13 8EX

01543 414704  
Sales@geneflow.co.uk  
Geneflow.co.uk