



***TransIT*[®]-COS Transfection Kit**

A high efficiency, low toxicity, DNA transfection reagent optimized for COS cell types

- **Unique Formulation**—Each kit contains a base and booster reagent to maximize transfection performance in COS cell types.
- **Low Cellular Toxicity**—Maintain cell density and reduce experimental biases.
- **High Efficiency Delivery**—Achieve 90-95% transfection efficiency in COS cells to ensure experimental success.

Data

Figure 1. The *TransIT*-COS Kit Achieves High Transfection Efficiency.

COS-7 cells were transfected with pEGFP using the *TransIT*-COS Kit in complete media for 24 hours.

