



***TransIT*®-Jurkat Transfection Reagent**

A high efficiency DNA transfection reagent optimized for Jurkat and other hard-to-transfect cell types

- **Unique Formulation**—Maximize transfection performance in Jurkat cells and other hard-to-transfect cells such as K562, RAW 264.7, THP-1 cells.
- **High Efficiency Delivery**—Achieve 25% transfection efficiency in Jurkat cells to ensure experimental success.
- **Deliver Single or Multiple Plasmids**—Suitable for many applications such as gene expression and promoter analysis.

Data

Figure 1. The *TransIT*-Jurkat Reagent Achieves High Transfection Efficiency.

Jurkat cells transfected with pEGFP using *TransIT*-Jurkat Reagent in complete media for 24 hours.

